



# InfinityLab LC/MSD Series Agilent LC/MSD XT

For OpenLAB CDS and OpenLAB CDS ChemStation Edition

## Data Sheet



The Agilent LC/MSD XT (liquid chromatograph/mass selective detector) achieves sensitivity and resolution specifications with autotune.

Parameter	Measure	Specification
Mass range		2–3,000 m/z
SIM sensitivity, ESI positive	1 pg of reserpine injected on column, quantifying the m/z 609 ion	Signal-to-noise ratio, RMS >300:1
SIM sensitivity, APCI positive	1 pg of reserpine injected on column, quantifying the m/z 609 ion	Signal-to-noise ratio, RMS >100:1
Maximum scan speed	Da/sec	10,400 u/sec (Da/sec)
Polarity switching	Switching time from positive to negative polarity	30 msec
Mass resolution (Autotune)	Full width at half maximum	Unit resolution ( $\leq 0.7$ u (Da))
Mass stability (drift)	Over 12 hours	$\leq 0.1$ u (Da) or 100 ppm, whichever is larger
Mass accuracy		$\pm 0.005$ u (Da) <sup>1</sup>
Spectral accuracy		99.0 % spectral accuracy <sup>1</sup>
Minimum SIM dwell time		5 msec
Maximum SIM transitions		100 SIM ions per time segment
Dynamic range		$> 6 \times 10^6$ (electronic)

<sup>1</sup> Only applicable with optional Cerno MassWorks Accurate Mass Software for LC/MSD using OpenLab CDS ChemStation Edition. Not verified during installation.  
u: Unified atomic mass unit



**Agilent Technologies**

## General system specifications

Parameter	Specification
Single point of control	Single-point data system method capability with full control of Agilent LCs including the 1100 Series, 1200 Series, Infinity Series and InfinityLab Series with the Agilent InfinityLab LC/MSD Series single quadrupole based systems.
Supported data systems	OpenLAB CDS, OpenLAB CDS ChemStation Edition, and MassHunter Data Analysis <sup>1</sup>
Time programming	<ul style="list-style-type: none"><li>• SIM ion selection</li><li>• Scan functions</li><li>• Fragmentor voltage</li><li>• External CAN valves</li><li>• MS diverter valve<sup>2</sup></li><li>• Drying gas flow</li><li>• Nebulizer pressure</li><li>• Drying gas temperature</li><li>• Capillary voltage (Positive)</li><li>• Capillary voltage (Negative)</li></ul>
Wide range of ionization sources	<ul style="list-style-type: none"><li>• Electrospray (ESI)</li><li>• Multimode source (simultaneous ESI and APCI)</li><li>• Atmospheric pressure chemical ionization (APCI)</li><li>• Atmospheric pressure photoionization (APPI)</li><li>• Agilent Jet Stream (AJS)</li></ul>
Detector	High energy conversion dynode and electron multiplier horn
Vacuum system	One turbo-molecular pump with one mechanical pump

<sup>1</sup> With data translation

<sup>2</sup> OpenLAB CDS only

## Ordering information

### G6135BA: Agilent LC/MSD XT System with OpenLAB CDS ChemStation Edition

Includes Agilent LC/MSD XT, based on single quadrupole technology, Electrospray ion source, and Agilent OpenLab ChemStation Edition data system with PC and display. System installation and familiarization included.

### G6135CA: Agilent LC/MSD XT System for OpenLAB CDS

Includes Agilent LC/MSD XT, based on single quadrupole technology, Electrospray ion source, and Agilent OpenLab CDS data system with PC and display. System installation and familiarization included.

### G2735L: Agilent Jet Stream Upgrade for the Agilent LC/MSD XT

### G6887AA: Cerno MassWorks Accurate Mass Software for LC/MSD

The above are bid specifications and are not performed at installation. See Site Preparation Guide, Installation Manual and Service Notes for additional product and specification information.

[www.agilent.com/chem/lcmsd](http://www.agilent.com/chem/lcmsd)

Agilent Technologies shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2017  
Published in the U.S.A.  
May 4, 2017  
5991-7830EN



**Agilent Technologies**