



Agilent 1260 Infinity Isocratic Pump

Features, Specifications and Ordering Details



HPLC workhorse for demanding QA/QC and routine applications

The optimized design of the Agilent 1260 Infinity Isocratic Pump ensures virtually pulse-free and stable solvent flow, with dual floating pistons in series, precisely servo-controlled. This delivery principle allows variable stroke volume optimized to selected flow rate. The broadened pressure range with maximum pressure of 600 bar allows the use of small particle size columns, longer columns, or alternative solvents of higher viscosity. It is the ideal pump for GPC/SEC applications where run-to-run and day-to-day precision of retention times is important.

Features

- 600 bar as new standard for additional usage of modern LC columns to obtain higher resolution and increased throughput.
- Low operating and maintenance costs – due to long instrument uptime.
- Future proof high investment protection since upgradeable to quaternary pump.
- Versatile wide flow range – from 0.05 mL/min up to 10 mL/min for analytical and semipreparative separations.
- Easy front access to all maintenance parts – parts can be exchanged quickly without the need to disassemble the instrument.
- Robust materials withstand most demanding applications – stainless steel, titanium, gold, ruby, sapphire, ceramics, PEEK and PTFE ensure long lifetime. Seal wash option for high salt mobile phases keeps maintenance to a minimum for lower operation costs.



Agilent Technologies

Specifications – Agilent 1260 Infinity Isocratic Pump

| Specifications Agilent 1260 Infinity Isocratic Pump (G1310B) | |
|--|---|
| Hydraulic system | Dual piston in series pump with servo-controlled variable stroke drive, floating pistons. |
| Settable flow range | Set points 0.001 - 5 mL/min, in 0.001 mL/min increments. |
| Flow range | 0.2 - 10.0 mL/min |
| Flow precision | ≤ 0.07 % RSD or ≤ 0.02 min SD, whatever is greater; based on retention time at constant room temperature. |
| Flow accuracy | ± 1 % or 10 µL/min whatever is greater, pumping degassed H ₂ O at 10 MPa (100 bar). |
| Pressure operating range | Operating range up to 60 MPa (600 bar, 8700 psi) up to 5 mL/min Operating range up to 20 MPa (200 bar, 2950 psi) up to 10 mL/min |
| Pressure pulsation | < 2 % amplitude (typically < 1.3 %) or < 0.3 MPa (3 bar), whatever is greater, at 1 mL/min isopropanol, at all pressures > 1 MPa (10 bar, 147 psi). |
| Compressibility compensation | User selectable, based on mobile phase compressibility. |
| Control | Agilent control software (e.g. ChemStation, EZChrom, OL, MassHunter). |
| Local control | Agilent Instant Pilot |
| Communications | Controller-area network (CAN), RS-232C, APG Remote: ready, start, stop and shut-down signals, LAN optional. |
| Safety and maintenance | Extensive diagnostics, error detection and display through Agilent LabAdvisor, leak detection, safe leak handling, leak output signal for shutdown of the pumping system. Low voltage in major maintenance areas. |
| GLP features | Early maintenance feedback (EMF) for continuous tracking of instrument usage in terms of seal wear and volume of pumped mobile phase with pre-defined and user settable limits and feedback messages. Electronic records of maintenance and errors. |
| Housing | All materials are recyclable. |

Ordering Details – Agilent 1260 Infinity Isocratic Pump

| Description | Product Number |
|---|----------------|
| Agilent 1260 Infinity Isocratic Pump. Maximum pressure 600 bar | G1310B |
| Tool kit for 1260/1290 Infinity LC | G1310B#001 |
| HPLC Starter-Kit incl. 0.17 mm ID capillaries | G1310B#002 |
| HPLC Starter-Kit incl. 0.12 mm ID capillaries | G1310B#003 |
| Lab Advisor | G1310B#004 |
| Active seal wash | G1398A |
| LAN Interface | G1310B#500 |
| Delete option for solvent-cabinet | G1310B#960 |
| Upgrade isocratic pump to quaternary | G4207A |
| Agilent 1260 Infinity Standard Degasser (<i>recommended</i>) | G1322A |
| Agilent 1260 Infinity Micro Degasser (<i>recommended</i>) | G1379B |

www.agilent.com/chem/1200

© Agilent Technologies, Inc., 2010
Published in USA, July 1, 2010
Publication Number 5990-6099EN



Agilent Technologies