


# Agilent 1260 Infinity Manual Injector

## Features and Ordering Details



### For low bench space footprint

Convenient modular accessory, dimensioned to match the footprint of all Agilent 1260 Infinity modules. Stack to create a compact, affordable manual injection liquid chromatography. It is compatible with up to 600-bar operating pressure.

### Features

- Lowest bench space requirements and highest flexibility with external HPLC columns by using the integrated column and base
- Long valve life, up to 600-bar operating pressure
- High instrument uptime through proven parts, e.g. ceramic rotor faces



**Agilent Technologies**

## Ordering Details – Agilent 1260 Infinity Manual Injector

Description	Product Number
Agilent 1260 Infinity manual injector with 20 $\mu$ L sample loop and mounting stand	G1328C
5- $\mu$ L sample loop	0101-1248
10- $\mu$ L sample loop	0100-1923
20- $\mu$ L sample loop	0100-1922
50- $\mu$ L sample loop	0100-1924
100- $\mu$ L sample loop	0100-1921
200- $\mu$ L sample loop	0101-1247
500- $\mu$ L sample loop	0101-1246
1-mL sample loop	0101-1245
2-mL sample loop	0101-1244
5-mL sample loop	0101-1243

[www.agilent.com/chem/1200](http://www.agilent.com/chem/1200)

© Agilent Technologies, Inc., 2010-2012

Published in USA, November 1, 2012

Publication Number 5990-6122EN



**Agilent Technologies**